

Pooled sunflower DNA

pAD2 - contains single unmethylated Pstl site.

Sample A	Buffer H	Invitrogen	15U/uL	2.7 uL	For 40 U
Sample B	Buffer O	Thermoscientific	10U/uL	4 uL	For 40 U
Sample C	Cutsmart	NEB	20U/uL	2	For 40U
Sample D	Cutsmart	NEB-HF	20U/uL	2	For 40U

Sample E	Buffer H	Invitrogen	15U/uL	2.7 uL	For 40 U
Sample F	Buffer O	Thermoscientific	10U/uL	4 uL	For 40 U
Sample G	Cutsmart	NEB	20U/uL	2	For 40U
Sample H	Cutsmart	NEB-HF	20U/uL	2	For 40U

Electrophoresis after 70 minutes - 40 uL of 200 uL loaded for each sample. Top lanes - A -D + uncut DNA. Bottom lanes - pAD2

A B C D Uncut

Lowest Is 500 to intense

GeneRule


E F G H

E F G H

Ladder Fermen

Lowest intense ladder marker Is 500 bp. Second lowest intense marker is 1500 bp.

GeneRuler 1 kb Plus DNA Ladder




Ladder = SM1331 -Fermentas/Thermofisher Electrophoresis after 18 hour digest - 60 uL of 160 uL loaded for each sample. Top lanes - A -D + uncut DNA. Bottom lanes - pAD2

D

G H

Uncut



Results - incomplete digestion after 70 minutes. 70 minute digestion is slightly greater than GBS protocol, prescribed, 3hrs with 10 U Pstl. Partial digestion observed overnight. None of the enzymes exhibit appreciable start activity at 40 U enzyme for 200 uL digestion.

Reaction setup from October 7, 2013

Testing Pstl NEB Pstl HF NEB Pstl Thermoscientific/Fermentas Pstl Invitrogen

Pool DNA from Dan B and Kate O DB209

DB183 DB238 JA52 DB320

Measured pooled DNA with nanodrop - 260/280; 260/230/ ng/uL

1.66 25.67 49.7

G5D MON014

Measured pooledDNA with nanodrop - 260/280; 260/230/ ng/uL

1.69 1.30 83.3

Digested 50 uL of DNA and 25 uL of DNA

Used volumes shown in column 4 to make 40 U of enzyme for each digest

Sample A	Buffer H	Invitrogen	15U/uL	2.7 uL	For 40 U
Sample B	Buffer O	Thermoscientific	10U/uL	4 uL	For 40 U
Sample C	Cutsmart	NEB	20U/uL	2	For 40U
Sample D	Cutsmart	NEB-HF	20U/uL	2	For 40U

Lot 0411106 Lot 0011010

200 uL digest 127.3 water Sample A 126 water B 128 water C 128 water D 20 buffer (see column 2 above) Pstl (see column 4 above) 2.7 A

4.0 B

2 C 2 D

DNA Samples A and B received 50 uL, C and D received 25 uL each