

## Sanger sequencing beads purification

### **SANGER SEQUENCING BEADS PURIFICATION**

(Basically identical to Elkin *et al.*, Biotechniques 2002)

This assumes that you start with a 10  $\mu$ l sequencing reaction in 0.2 ml PCR tubes.

- add 20  $\mu$ l of BET solution, mix well by pipetting
- Incubate at room temperature for 10 minutes
- Place tubes a plate magnet for about 1 minute
- Remove the supernatant
- Add 30  $\mu$ l of 70% ethanol
- Remove the ethanol (be careful to remove all of it)
- Air dry for 5-10 minutes
- Remove from the magnet and add 15  $\mu$ l of milliQ water
- Incubate at room temperature for 5 minutes
- Place tubes on a plate magnet for 2 minutes
- Transfer 14  $\mu$ l of water solution to new tubes

#### **BET solution (2 ml, 100 cleanups)**

- |                         |             |
|-------------------------|-------------|
| - Ethanol 100%          | 1.6 ml      |
| - milliQ water          | 190 $\mu$ l |
| - Tetra Ethylene Glycol | 160 $\mu$ l |
| - Sera-Mag SpeedBeads   | 50 $\mu$ l  |

The original protocol does not say anything about washing the Sera-Mag beads before adding them to the BET solution. For good measure (to remove the sodium azide), I put the 50  $\mu$ l on a magnet, remove the supernatant, was once with Tris 10 mM pH 8.0, and then re-suspend them in 50  $\mu$ l of Tris 10 mM pH 8.0